

02

CUVETTE & SAMPLE CUP

| Classification | REF. | Page |
|----------------|------|-------|
| ◎ CUVETTE | 081 | 17-20 |
| ◎ SAMPLE CUP | 082 | 21-24 |

- ◎ KANG JIAN® cuvette & Sample Cup is made of imported high polymer PS material.
- ◎ High precision mould makes the cuvettes smooth, bright and clean.
- ◎ Accurate size and structure makes the cuvette adapted to different coagulometer, biochemical analyzer, and multiple analytical systems.

KANGJIAN Medical Apparatus CO.
NO.16 Zhangqian Road • Jiangyan Taizhou
• Jiangsu Province • China

Customer service, sales and technical support
Phone: +86 523 82282698 • +86 523 88857012
Fax: +86 523 88665588 • +86 523 88857712
E-mail: export@china-kangjian.com

Cuvette



081-101
(KJ133)



081-102
(KJ134)



081-103
(KJ135)



081-104
(KJ136)



081-105
(KJ139)



081-106
(KJ111)



081-107
(KJ111-1)



081-108
(KJ111-2)



081-109
(KJ112-1)



081-110
(KJ112-2)



081-111
(KJ104)



081-112
(KJ104-1)

| REF | Original REF. | MAT | Size | Box | Tray | Info |
|---------|---------------|-----|----------------|------|-------|-------------------------|
| 081-101 | KJ133 | PS | 12.5×12.5×45mm | 100 | 2000 | Application: 340-800mm |
| 081-102 | KJ134 | PS | 12.5×12.5×45mm | 100 | 2000 | Application: 340-800mm |
| 081-103 | KJ135 | PS | 12.5×12.5×45mm | 100 | 2000 | Application: 340-800mm |
| 081-104 | KJ136 | PS | 12.5×12.5×45mm | 100 | 2000 | Application: 340-800mm |
| 081-105 | KJ139 | PS | 15.7×11.8×46mm | 100 | 2000 | |
| 081-106 | KJ111 | PS | ∅15×35mm | 2000 | 20000 | For Sysmex CA 530 |
| 081-107 | KJ111-1 | PS | ∅15×35mm | 1000 | 24000 | For Sysmex CA 7000 |
| 081-108 | KJ111-2 | PS | ∅15×35mm | 3000 | 24000 | For Sysmex CA 1500 |
| 081-109 | KJ112-1 | PS | ∅15×35mm | 2000 | 24000 | For Sysmex CA 2000i |
| 081-110 | KJ112-2 | PS | ∅8.8×30mm | 2000 | 20000 | For SEKISUI |
| 081-111 | KJ104 | PS | 60×32mm | 29 | 2001 | For German BE Automatic |
| 081-112 | KJ104-1 | PS | 119×32mm | 8 | 1160 | For German BE Automatic |

Cuvette



| REF | Original REF. | MAT | REF | | | | |
|---------|---------------|-----|-------------|------|--|-------|----------------------------|
| 081-113 | KJ101 | PS | 7,5×8×23mm | 1000 | | 25000 | For German BE 1 Channel |
| 081-114 | KJ102 | PS | 7,5×29×21mm | 600 | | 6000 | For German BE 2 Channel |
| 081-115 | KJ103 | PS | 7,5×61×21mm | 300 | | 6000 | For German BE 4 Channel |
| 081-116 | KJ105 | PS | 11×30×27mm | 100 | | 4000 | For German BECL4 2 Channel |
| 081-117 | KJ106 | PS | 11×60×27mm | 50 | | 1000 | For German BECL4 4 Channel |
| 081-118 | KJ107 | PS | 11×120×27mm | 25 | | 500 | For German BECL4 8 Channel |
| 081-119 | KJ108 | PS | Ø17×26mm | 1000 | | 10000 | For German TECO |
| 081-120 | KJ109 | PS | 11×31×19mm | 1000 | | 10000 | For German TECO 2 Channel |
| 081-121 | KJ110 | PS | Ø10×19mm | 1000 | | 25000 | For German TECO 1 Channel |
| 081-122 | KJ110-1 | PS | Ø10×23mm | 1000 | | 25000 | For German TECO 1 Channel |
| 081-123 | KJ114 | PS | 7,5×8×23mm | 1000 | | 10000 | For Spain Cruor Apparatus |
| 081-124 | KJ115 | PS | 11×120×27mm | 1000 | | 10000 | For German LABOR |

Cuvette



081-125
(KJ124)



081-126
(KJ125)



081-127
(KJ132)



081-128
(KJ137)




081-129
(KJ138)



081-130
(KJ116)



081-131
(KJ118)



081-132
(KJ122)



081-133
(KJ120)



081-134
(KJ121)


| REF | Original REF. | MAT | REF | REF | REF | REF | REF |
|---------|---------------|-----|--------------|------|-----|-------|-------------------------------|
| 081-125 | KJ124 | PS | 11×60×27mm | 50 | | 3000 | For Beijing Perlong |
| 081-126 | KJ125 | PS | 11×30×27mm | 25 | | 500 | For Swiss SEROZYME |
| 081-127 | KJ132 | PS | Ø11×23mm | 1000 | | 10000 | For Sigma KC-40 |
| 081-128 | KJ137 | PS | Ø12×30mm | 1000 | | 10000 | FC Serum Cup |
| 081-129 | KJ138 | PS | Ø10×30mm | 1000 | | 10000 | FC Serum Cup |
| 081-130 | KJ116 | PS | 47,5×9×27mm | 200 | | 4000 | For American Pacific |
| 081-131 | KJ118 | PS | 45×19×28mm | 100 | | 2000 | For Zhongshan Peikang |
| 081-132 | KJ122 | PS | Ø108mm | 1 | | 500 | For American Coulter |
| 081-133 | KJ120 | PS | 66,5×16×24mm | 175 | | 2100 | For Precil / For French STAGO |
| 081-134 | KJ121 | PS | 66×16×30mm | 140 | | 2800 | For Steellex |

Cuvette


081-135
(KJ123)



081-136
(KJ127)



081-137
(KJ128)



081-138
(KJ129)



081-139
(KJ130)



REF

Original REF.

MAT



081-135

KJ123

PS

96×9×18mm

2000

For Beckman

081-136

KJ127

PS

96×9×18mm

400

For Roche

081-137

KJ128

PS

1200

For Specific Protein Analyzer

081-138

KJ129

PS

2000

For BN100 Specific Protein Analyzer

081-139

KJ130

PS

1200

For French Biomerieux

Sample Cup



082-101
(KJ206)



082-102
(KJ211)



082-103
(KJ212)



082-104
(KJ213)



082-105
(KJ214)



082-106
(KJ215)



082-107
(KJ216)



082-108
(KJ219)



082-109
(KJ220)



082-110
(KJ224)



082-111
(KJ225)



082-112
(KJ226)



082-123
(KJ241)








082-124
(KJ242)

| REF | Original REF. | MAT | UNITS | PKT | BOX | TRAY | INFO |
|---------|---------------|---------|------------|------|-------|------|--|
| 082-101 | KJ206 | PS + PP | Ø24×45mm | 180 | 1800 | | For Swiss AC900 series |
| 082-102 | KJ211 | PS | Ø16×38mm | 500 | 5000 | | For Hitachi 7150, 7060 series |
| 082-103 | KJ212 | PS | Ø13×24mm | 1000 | 10000 | | For American Beckman CX series |
| 082-104 | KJ213 | PS | Ø14×25mm | 1000 | 10000 | | For American Beckman 700/American Dupont |
| 082-105 | KJ214 | PS | Ø13×20mm | 1500 | 15000 | | For Kodak |
| 082-106 | KJ215 | PS | Ø13×17mm | 1500 | 13500 | | For American IL |
| 082-107 | KJ216 | PS | Ø12×20mm | 2000 | 20000 | | For American Corning 550, 560 series |
| 082-108 | KJ219 | PS | Ø15×40mm | 500 | 5000 | | For American Abbott CCX |
| 082-109 | KJ220 | PS | Ø13×38mm | 500 | 6000 | | For Hitachi Minimum |
| 082-110 | KJ224 | PS | | 1000 | 9000 | | For Hitachi (Conical sample cup) |
| 082-111 | KJ225 | PS | Ø10×49mm | 500 | 5000 | | For American Abbott |
| 082-112 | KJ226 | PS | Ø17.5×53mm | 400 | 4000 | | For American Abbott AXSYM |
| 082-123 | KJ241 | PS | | | 2000 | | DiaSorin biochemical Cup (six hole) |
| 082-124 | KJ242 | PS | | | 1000 | | For Dade Behring protein analyzer |

Sample Cup



| REF | Original REF. | MAT |  |  |  |  |  |
|---------|---------------|-----|---|---|---|---|---|
| 082-113 | KJ221 | PS | Ø7×26mm | 2000 | | 20000 | For Bowlinman / Roche |
| 082-114 | KJ222 | PP | Ø9×18mm | 2000 | | 20000 | For Shanghai AnTai |
| 082-115 | KJ223 | PP | Ø11×19mm | 2000 | | 20000 | For American M10 |
| 082-116 | KJ223-1 | PP | Ø15×22mm | 1000 | | 10000 | For SINNOWA |
| 082-117 | KJ246 | PP | Ø15×25mm | 1000 | | 10000 | For Simens |
| 082-118 | KJ236 | PP | 40×11mm | 1000 | | 10000 | For Abbott Architect i2000 SR |
| 082-119 | KJ236-1 | PP | | 1000 | | 10000 | For Autobio |
| 082-120 | KJ238 | PP | | 1000 | | 10000 | For Beckman DXI800 |
| 082-121 | KJ239 | PP | | 1000 | | 10000 | For Simens |
| 082-122 | KJ247 | PS | | 1000 | | 10000 | For Sysmex |
| 082-129 | KJ217 | PP | 8×34mm | 1000 | | 10000 | For Roche COBAS |
| 082-125 | KJ227 | PS | 37×14×21mm | 200 | | 6000 | For Ciba Corning 550 Express |
| 082-126 | KJ228 | PS | 53×9×18mm | 25 | | 2500 | For Finland Cornyi |

Sample Cup



082-127
(KJ229)



082-128
(KJ230)



Type A

A-Type Cap
Quantitative Suction Dropper



Type B

B-Type Cap
Quantitative Suction Dropper



Type C

C-Type Cap
Quantitative Suction Dropper



082-137
(KJ240)








082-138
(KJ240-1)





082-139
(KJ240-2)








082-140
(KJ240-3)

| REF | Original REF. | MAT |  |  |  |  |  |
|---------|---------------|-----|---|---|---|---|---|
| 082-127 | KJ229 | PS | | 15 | | 1050 | For American MD2000 / 1800 |
| 082-128 | KJ230 | PS | | 15 | | 1050 | For American MD3000 |
| 082-137 | KJ240 | PS | 5µl | 1000 | | 10000 | Quantitative Suction Dropper, For Type A / B / C |
| 082-138 | KJ240-1 | PS | 10µl | 1000 | | 10000 | Quantitative Suction Dropper, For Type A / B / C |
| 082-139 | KJ240-2 | PS | 7.5µl | 1000 | | 10000 | Quantitative Suction Dropper, For Type A / B / C |
| 082-140 | KJ240-3 | PS | 10µl | 1000 | | 10000 | Quantitative Suction Dropper, For Type A / B / C |

Sample Cup



| REF | Original REF. | MAT |  |  |  |  |  |
|---------|---------------|-----|---|---|---|---|---|
| 082-131 | KJ231 | | 138×121×54mm | 48 | 96 | Set Packed Sample Cup & Tip for Roche | |
| 082-132 | KJ232 | | 112×94×45mm | 30 | 120 | Set Packed Tip for Roche | |
| 082-133 | KJ233 | | 95×59×30mm | 60 | 240 | Set Packed Sample Cup for Roche | |
| 082-136 | KJ245 | PS | 18×9×7mm | 100 | 2000 | For Beckman ACL-TOP | |

